

Circular Sample Cutter

manual

Code: D.513.xxx



Usage

To prepare round samples for the grammage determination.

Features

Sharp, circularly guided knives achieve true-to-size round samples with absolutely smooth edges in the drawing cut. Even critical materials such as thin foils, tissue paper, corrugated cardboard, hosiery fabrics, nonwovens, carpets, artificial leather, fluffy or velvety fabrics can be cut without any issues. The cutting base is a poro rubber plate, into which the blades can penetrate without danger. The standard device has a cutting depth of 5 mm, but some models are also offered with a cutting depth of 10 mm.

Delivery content

- Round sample cutter
- Poro rubber base plate
- User manual
- Certificates

Models

Code	Sample size
D.513.010	10 cm ² = Ø 35,7 mm
D.513.025	25 cm ² = Ø 56,4 mm
D.513.050	50 cm ² = Ø 79,8 mm
D.513.100	100 cm ² = Ø 112,8 mm
D.513.D47	17.34cm ² = Ø 47 mm
D.513.D50	19.63 cm ² = Ø 50 mm

Customized sizes are available!
 For more information, contact your personal distributor or PTE Austria directly!